

BIOTHANE™ G1 STRAPS

BIOTHANE® G1 STRETCHER/COT STRAPS



BIOTHANE® G1 backboard straps are made to secure patients to backboard.

For a custom fit, economy disposable poly straps are easily adjustable.

2" **BIOTHANE® G1** webbing (1400 pound break strength)

Meets FDA requirements.



INSTRUCTIONS:

- 1.) Inspect strap before use.
- 2.) Position patient.
- 3.) Strap should have proper tension with no twists.
- 4.) Patient should be flat on backboard/stretchers.
- 5.) Wrap strap around patient and backboard/stretchers.
- 6.) Buckle and tighten strap to secure.

Intended as a guide only.
 Consult your local medical agency for all protocols.

1B Metal Push Button Buckle 5000 Pound Break Strength	7B Plastic Side Release Buckle 400 Pound Break Strength
2E Metal Swivel Speed Clip Ends 500 Pound Break Strength	6E Metal Loop Ends 495 Pound Break Strength

Black	Orange	Neon Yellow	Blue
	Silver	Golden Yellow	



For use by qualified and trained personnel only.