

## HEAD IMMOBILIZERS

### HEAD GUARD™



**Item # 50100**

Head Guard™ Head Immobilizer, each

**Item # 50106**

Head Guard™ Head Immobilizer, case/6



The **Head Guard™** is a reusable custom tailored head immobilizer that fits most spine boards. The side blocks (held in place with extra strong Hook & Loop that works even in water) are adjustable, shifting to accommodate almost any angle. Large ear holes make it easy for patient evaluation. The head and chin straps hold the head securely in place. Vinyl dipped foam blocks provide firm but resilient support and a clean hygienic surface. Easily attaches to spine boards with all Hook & Loop strapping.

#### DIMENSIONS:

- 15 ½" X 10" Base
- 6" H X 9½" W Tapered Pillows
- 39" Head & Chin Straps
- 25" attachment straps

#### AVAILABLE COLORS:

<b>Blue</b>	<b>Black</b>
<b>Orange</b>	<b>Red</b>



**Head Guard™  
Pillows, Base and Head & Chin Straps  
ALSO SOLD SEPERATELY!**

Item #50101


Head Guard™ Head Immobilizer Pillows, pair

Item #50102

Head Guard™ Head Immobilizer Base, each

Item #50103

Head Guard™ Head Immobilizer Head & Chin Straps, pair



**Item #50104BL**

Pediatric Pillows Also Available

