

HEAD IMMOBILIZERS

HEAD STAY™



Item # 50000

The **Head Stay™** is an easy way to stabilize head or neck injuries. The Head Stay™ is made of water resistant cardboard. One blue adhesive head strap included. Each packaged flat for easy storage.

- 15" X 5" Base
- 15" head stay panels - Both sides of immobilizer are 6" when folded
- Can be easily adjusted to fit most patients
- Blue center foam pad for head comfort
- Highly visible **red** sticky pull tabs that secure firmly to board

