

SOFT STRETCHERS

SCOOP UP™



Item #30070YL

Scoop Up™ Soft Stretcher, Yellow/Orange, each

Item #30070BK

Scoop Up™ Soft Stretcher, Black/Orange, each

The **Scoop Up™ Stretcher** is designed to transport patients from confined areas. The patient is carried in an upright sitting position with 6 handholds. It measures 20" X 34" with 2" continuous and reinforced yellow or black nylon straps. The Scoop Up™ can be easily stored.

Made with Vinyl for easy clean and disinfecting.

