

The background of the entire image is a close-up, high-resolution photograph of the United States flag. The flag is shown waving, with the blue field of stars in the upper left and the red and white stripes filling the rest of the frame. The lighting creates a sense of depth and texture in the fabric.

DMS

A red ECG (heart rate) line graphic that starts under the 'D' of 'DMS', dips, and then rises to underline the 'S'.

Dick Medical Supply

TABLE OF CONTENTS

TERMS AND CONDITIONS	2
NYLON STRETCHER/COT STRAPS.....	3&4
POLYPROPYLENE STRETCHER/COT STRAPS.....	5&6
IMPERVIOUS VINYL STRETCHER/COT STRAPS.....	7&8
CROSS HARNESS SYSTEMS.....	9
SHOULDER RESTRAINT HARNESS SYSTEMS	10
BIOTHANE® G1.....	11
DISPOSABLE POLYPROPYLENE STRAPS	12
ECONOMY DISPOSABLE RED, WHITE & BLUE STRAP SETS	12
BODY STRAP SYSTEMS.....	13&14
DMS BRAND KITS	15&16
PPE	17&18
CAPS AND COVERS	18
BACKBOARDS	19
PEDIATRIC CHILD RESTRAINT.....	19
CARRY POUCH™	19
SOFT STRETCHERS	20
HIP & PELVIC STABILIZERS	20
GAIT BELTS	21&22
LIMB RESTRAINTS	23
STAIR CHAIR REPLACEMENT STRAPS.....	24
EQUIPMENT HOLDING & CINCH STRAPS	24
LEND-A-HAND STRETCHER GRIPS	24
HANDS OFF™ SECURITY BUCKLE	24
HEAD IMMOBILIZERS.....	25&26
HARE TRACTION SPLINT REPLACEMENT STRAPS.....	26
PADDED BOARD SPLINTS.....	27
CARDBOARD SPLINTS	27
CARCINOGEN CONTAINMENT GEAR BAG	28
CARGO NET	28
IV ARMBOARDS.....	28
PREMIUM & ECONOMY COTS	29
RESCUE CHAIR.....	29
DMS.....	30

CUSTOM SEWING & MANUFACTURING

Our extensive industrial sewing abilities allow us to produce items for a wide range of industries. These include military, industrial, medical, sports and many more. For any custom sewing applications or custom product needs contact DMS today @ 614-444-2300.

DAVID DICK - ddick@dms4ems.com

GREG KNEELAND - gkneeland@dms4ems.com

KATRINA COPELAND - kcopeland@dms4ems.com

TERMS AND CONDITIONS

Business Hours:

Our customer service staff is ready to assist you with your order or any questions you may have Monday through Friday from 7:00 AM to 3:30 PM Eastern Time. Contact (614) 444-2300 or by Email. David Dick - ddick@dms4ems.com, Greg Kneeland - gkneeland@dms4ems.com or Katrina Copeland - kcopeland@dms4ems.com

Terms:

Net 30 Days for all established accounts.

Payments:

On all Purchase Orders submitted, purchaser has the responsibility for proper and timely payments. Any account over 30 days will be charged 2% per month (24% annual). Returned checks are subject to a \$35.00 process charge and a \$20.00 rebilling charge. DMS does accept credit cards. (MasterCard, Visa and American Express) A credit card fee will be charged for ALL credit card transactions.

Pricing:

All prices are subject to change without written and/or verbal notice.

Return Policy:

- Returns must be approved by Dick Medical Supply with a return authorization number.
- No cash refunds. Credit on account will be issued.
- There is a 90 day return policy on all cardboard splints.
- Any product that is customized, custom labeled, screen printed or custom ordered **can not** be returned.
- All products returned must be in original package.
- All products will be inspected upon return. If any damages are found during the inspection, DMS will contact you with the description of the damage and the quantity. Any damaged product will not be accepted and therefore will not be issued credit.
- A 20% restocking fee will apply for all returned products.

Freight and Shipping Policies:

All orders submitted to Dick Medical are shipped FOB Columbus, Ohio, unless otherwise specified by the purchaser. Dick Medical's current carriers are UPS and FedEx, with all shipping charges prepaid and then added to the purchaser's invoice. If you have a UPS and FedEx Account number, please make a note on your purchase order. Any shipments that are shipped on customer accounts will be required to file a claim with such carrier for any damaged/lost cartons or items.

Third Party & Customer Designated Shippers:

If you provide a 3rd party account number or designate a specific carrier, you assume responsibility from the time the order is ready for pickup and are responsible for any lost or damaged goods.

Warranties:

No Express warranties and no implied warranties on products sold.

Satisfaction Guarantee:

Our goal at Dick Medical Supply is to supply you with high quality products, great customer service and timely shipping. Should a problem arise, we will make every effort to correct the problem in a timely manner. If you need any assistance, please contact us at (614) 444-2300.

Rush and/or Same Day Order Shipping:

\$50.00 fee will be applied.



DMS offers a large variety of stretcher and cot straps. 2" NYLON webbing - 4000 POUND BREAK STRENGTH
Contact us for further information.



1B METAL PUSH BUTTON	9B PLASTIC CAM
2B METAL CAM	1E LOOP END
3B METAL ROLLER	2E METAL SWIVEL SPEED CLIP
4B METAL DROP JAW	3E NON-METALSWIVEL SPEED CLIP
5B METAL DOUBLE D-RING	4E BIG BRASS SWIVEL SPEED CLIP
6B PLASTIC LADDERLOC	5E PLASTIC SWIVEL SPEED CLIP
7B PLASTIC SIDE RELEASE	6E METAL LOOP
8B	

Metal Push Button Buckle, 7', 1 piece, each

Metal Push Button Buckle, 9', 1 piece, each

Metal Push Button Buckle w/ Loop Ends, 5', 2 piece, each

Metal Push Button Buckle w/ Loop Ends, 7', 2 piece, each

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 5', 2 piece, each

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 7', 2 piece, each

Metal Push Button Buckle w/ Non-Swivel Speed Clip Ends, 5', 2 piece, each

Metal Push Button Buckle w/ Non-Swivel Speed Clip Ends, 7', 2 piece, each

Metal Push Button Buckle w/ Big Brass Swivel Clip Ends, 5', 2 piece, each

Metal Push Button Buckle w/ Big Brass Swivel Clip Ends, 7', 2 piece, each

Metal Push Button Buckle w/ Metal Loop Ends, 5', 2 piece, each

Metal Push Button Buckle w/ Metal Loop Ends, 7', 2 piece, each

Metal Push Button Buckle Extension Strap, 2', 1 piece, each

Metal Push Button Buckle Extension Strap, 3', 1 piece, each

Metal Push Button Buckle Extension Strap, 4', 1 piece, each

Metal Cam Buckle, 7', 1 piece, each

Metal Cam Buckle, 9', 1 piece, each

Metal Cam Buckle w/ Loop Ends, 5', 2 piece, each

Metal Cam Buckle w/ Loop Ends, 7', 2 piece, each



Hands Off™ Security Buckle available.
CALL FOR FURTHER INFORMATION!

BUCKLE SHIELD™ Security Buckle available.
CALL FOR FURTHER INFORMATION!



NYLON STRETCHER / COT STRAPS

METAL ROLLER BUCKLE

Item # 13071

Metal Roller Buckle, 7', 1 piece, each

Item # 13091

Metal Roller Buckle, 9', 1 piece, each

METAL DROP JAW BUCKLE

Item # 14091

Metal Drop Jaw Buckle, 9', 1 piece, each

Item # 14152

Metal Drop Jaw Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 14172

Metal Drop Jaw Buckle, w/ Loop Ends, 7', 2 piece, each

METAL DOUBLE D-RING

Item # 15091

Metal Double D-Ring Buckle, 9', 1 piece, each

Item # 150121

Metal Double D-Ring Buckle, 12', 1 piece, each

Item # 150151

Metal Double D-Ring Buckle, 15', 1 piece, each

PLASTIC LADDERLOC BUCKLE

Item # 16071

Plastic Ladderloc Buckle, 7', 1 piece, each

Item # 16091

Plastic Ladderloc Buckle, 9', 1 piece, each

PLASTIC SIDE RELEASE BUCKLE

Item # 17071

Plastic Side Release Buckle, 7', 1 piece, each

Item # 17091

Plastic Side Release Buckle, 9', 1 piece, each

Item # 17152

Plastic Side Release Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 17172

Plastic Side Release Buckle, w/ Loop Ends, 7', 2 piece, each

Item # 17252

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 17272

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 7', 2 piece, each

Item # 17352

Plastic Side Release Buckle, w/ Non-Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 17372

Plastic Side Release Buckle, w/ Non-Metal Swivel Speed Clip Ends, 7', 2 piece, each

Item # 17452

Plastic Side Release Buckle, w/ Big Brass Swivel Speed Clip Ends, 5', 2 piece, each

Item # 17472

Plastic Side Release Buckle, w/ Big Brass Swivel Speed Clip Ends, 7', 2 piece, each

PLASTIC DOUBLE ADJUST BUCKLE

Item # 18152

Plastic Double Adjust Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 18172

Plastic Double Adjust Buckle, w/ Loop Ends, 7', 2 piece, each

Item # 18252

Plastic Double Adjust Buckle, w/ Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 18272

Plastic Double Adjust Buckle, w/ Metal Swivel Speed Clip Ends, 7', 2 piece, each

PLASTIC CAM BUCKLE

Item # 19071

Plastic Cam Buckle, 7', 1 piece, each

Item # 19091

Plastic Cam Buckle, 9', 1 piece, each

AVAILABLE NYLON COLORS:

Black	Blue	Green
Maroon	Olive Drab	Orange
Red	Yellow	



Strap Carrying Case

The **Strap Carrying Case** is made with lined Heavy Duty Vinyl that is reinforced with 1" Nylon webbing. Velcro® closure. Belt loops on the back for easy carrying. Holds up to two straps.

Item #60004BK



SCREEN-PRINTING & EMBROIDERY AVAILABLE

Call for further information.

614-444-2300 | www.dickmedical.com

POLYPROPYLENE STRETCHER / COT STRAPS

DMS offers a large variety of stretcher and cot straps. 2" NYLON webbing - 4000 POUND BREAK STRENGTH
Contact us for further information.

Buckles & Ends



1B METAL PUSH BUTTON	9B PLASTIC CAM
2B METAL CAM	1E LOOP END
3B METAL ROLLER	2E METAL SWIVEL SPEED CLIP
4B METAL DROP JAW	3E NON-METAL SWIVEL SPEED CLIP
5B METAL DOUBLE D-RING	4E BIG BRASS SWIVEL SPEED CLIP
6B PLASTIC LADDERLOC	5E PLASTIC SWIVEL SPEED CLIP
7B PLASTIC SIDE RELEASE	6E METAL LOOP
8B PLASTIC DOUBLE AJUST	

METAL PUSH BUTTON BUCKLE

Item # 21071

Metal Push Button Buckle, 7', 1 piece, each

Item # 21091

Metal Push Button Buckle, 9', 1 piece, each

Item # 21152

Metal Push Button Buckle w/ Loop Ends, 5', 2 piece, each

Item # 21172

Metal Push Button Buckle w/ Loop Ends, 7', 2 piece, each

Item # 21252

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 5', 2 piece, each

Item # 21272

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 7', 2 piece, each

Item # 21352

Metal Push Button Buckle w/ Non-Swivel Speed Clip Ends, 5', 2 piece, each

Item # 21372

Metal Push Button Buckle w/ Non-Swivel Speed Clip Ends, 7', 2 piece, each

Item # 21452

Metal Push Button Buckle w/ Big Brass Swivel Clip Ends, 5', 2 piece, each

Item # 21472

Metal Push Button Buckle w/ Big Brass Swivel Clip Ends, 7', 2 piece, each

Item # 21652

Metal Push Button Buckle w/ Metal Loop Ends, 5', 2 piece, each

Item # 21672

Metal Push Button Buckle w/ Metal Loop Ends, 7', 2 piece, each

EXTENSION STRAPS

Item # 21021

Metal Push Button Buckle Extension Strap, 2', 1 piece, each

Item # 21031

Metal Push Button Buckle Extension Strap, 3', 1 piece, each

Item # 21041

Metal Push Button Buckle Extension Strap, 4', 1 piece, each

METAL CAM BUCKLE

Item # 22071

Metal Cam Buckle, 7', 1 piece, each

Item # 22091

Metal Cam Buckle, 9', 1 piece, each

Item # 22152

Metal Cam Buckle w/ Loop Ends, 5', 2 piece, each

Item # 22172

Metal Cam Buckle w/ Loop Ends, 7', 2 piece, each



Hands Off™ Security Buckle available.
CALL FOR FURTHER INFORMATION!

BUCKLE SHIELD™ Security Buckle available.
CALL FOR FURTHER INFORMATION!



POLYPROPYLENE STRETCHER / COT STRAPS

METAL ROLLER BUCKLE

Item # 23071

Metal Roller Buckle, 7', 1 piece, each

Item # 23091

Metal Roller Buckle, 9', 1 piece, each

METAL DROP JAW BUCKLE

Item # 24091

Metal Drop Jaw Buckle, 9', 1 piece, each

Item # 24152

Metal Drop Jaw Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 24172

Metal Drop Jaw Buckle, w/ Loop Ends, 7', 2 piece, each

METAL DOUBLE D-RING

Item # 25091

Metal Double D-Ring Buckle, 9', 1 piece, each

Item # 250121

Metal Double D-Ring Buckle, 12', 1 piece, each

Item # 250151

Metal Double D-Ring Buckle, 15', 1 piece, each

PLASTIC LADDERLOC BUCKLE

Item # 26071

Plastic Ladderloc Buckle, 7', 1 piece, each

Item # 26091

Plastic Ladderloc Buckle, 9', 1 piece, each

PLASTIC SIDE RELEASE BUCKLE

Item # 27071

Plastic Side Release Buckle, 7', 1 piece, each

Item # 27091

Plastic Side Release Buckle, 9', 1 piece, each

Item # 27152

Plastic Side Release Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 27172

Plastic Side Release Buckle, w/ Loop Ends, 7', 2 piece, each

Item # 27252

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 27272

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 7', 2 piece, each

Item # 27352

Plastic Side Release Buckle, w/ Non-Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 27372

Plastic Side Release Buckle, w/ Non-Metal Swivel Speed Clip Ends, 7', 2 piece, each

Item # 27452

Plastic Side Release Buckle, w/ Big Brass Swivel Speed Clip Ends, 5', 2 piece, each

Item # 27472

Plastic Side Release Buckle, w/ Big Brass Swivel Speed Clip Ends, 7', 2 piece, each

PLASTIC DOUBLE ADJUST BUCKLE

Item # 28152

Plastic Double Adjust Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 28172

Plastic Double Adjust Buckle, w/ Loop Ends, 7', 2 piece, each

Item # 28252

Plastic Double Adjust Buckle, w/ Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 28272

Plastic Double Adjust Buckle, w/ Metal Swivel Speed Clip Ends, 7', 2 piece, each

PLASTIC CAM BUCKLE

Item # 29071

Plastic Cam Buckle, 7', 1 piece, each

Item # 29091

Plastic Cam Buckle, 9', 1 piece, each

AVAILABLE POLYPROPYLENE COLORS:

Black	Blue	Green
Neon Green	Orange	Red
White	Yellow	

1" METAL SQUEEZE BUCKLE

Item # 24021

Black Polypropylene, 1" Metal Squeeze Buckle, 2', 1 piece, each

Item # 24031

Black Polypropylene, 1" Metal Squeeze Buckle, 3', 1 piece, each

Item # 24041

Black Polypropylene, 1" Metal Squeeze Buckle, 4', 1 piece, each

Item # 24061

Black Polypropylene, 1" Metal Squeeze Buckle, 6', 1 piece, each



Strap Carrying Case

The **Strap Carrying Case** is made with lined Heavy Duty Vinyl that is reinforced with 1" Nylon webbing. Velcro® closure. Belt loops on the back for easy carrying. Holds up to two straps.

Item #60004BK

Strap Carrying Case, each



SCREEN-PRINTING & EMBROIDERY AVAILABLE

Available in all webbing categories. Call for further information.
614-444-2300 |

IMPERVIOUS VINYL STRETCHER / COT STRAPS

DMS offers a large variety of stretcher and cot straps. 2" NYLON webbing - **4000 POUND BREAK STRENGTH**
Contact us for further information.

Buckles & Ends



1B METAL PUSH BUTTON	9B PLASTIC CAM
2B METAL CAM	1E LOOP END
3B METAL ROLLER	2E METAL SWIVEL SPEED CLIP
4B METAL DROP JAW	3E NON-METALSWIVEL SPEED CLIP
5B METAL DOUBLE D-RING	4E BIG BRASS SWIVEL SPEED CLIP
6B PLASTIC LADDERLOC	5E PLASTIC SWIVEL SPEED CLIP
7B PLASTIC SIDE RELEASE	6E METAL LOOP
8B PLASTIC DOUBLE AJUST	

METAL PUSH BUTTON BUCKLE

Item # 31071

Metal Push Button Buckle, 7', 1 piece, each

Item # 31091

Metal Push Button Buckle, 9', 1 piece, each

Item # 31152

Metal Push Button Buckle w/ Loop Ends, 5', 2 piece, each

Item # 31172

Metal Push Button Buckle w/ Loop Ends, 7', 2 piece, each

Item # 31252

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 5', 2 piece, each

Item # 31272

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 7', 2 piece, each

Item # 31352

Metal Push Button Buckle w/ Non-Swivel Speed Clip Ends, 5', 2 piece, each

Item # 31372

Metal Push Button Buckle w/ Non-Swivel Speed Clip Ends, 7', 2 piece, each

Item # 31452

Metal Push Button Buckle w/ Big Brass Swivel Clip Ends, 5', 2 piece, each

Item # 31472

Metal Push Button Buckle w/ Big Brass Swivel Clip Ends, 7', 2 piece, each

Item # 31652

Metal Push Button Buckle w/ Metal Loop Ends, 5', 2 piece, each

Item # 31672

Metal Push Button Buckle w/ Metal Loop Ends, 7', 2 piece, each

EXTENSION STRAPS

Item # 31021

Metal Push Button Buckle Extension Strap, 2', 1 piece, each

Item #31031

Metal Push Button Buckle Extension Strap, 3', 1 piece, each

Item # 31041

Metal Push Button Buckle Extension Strap, 4 , 1 piece, each

METAL CAM BUCKLE

Item # 32071

Metal Cam Buckle, 7', 1 piece, each

Item # 32091

Metal Cam Buckle, 9', 1 piece, each

Item # 32152

Metal Cam Buckle w/ Loop Ends, 5', 2 piece, each

Item # 32172

Metal Cam Buckle w/ Loop Ends, 7', 2 piece, each



Hands Off™ Security Buckle available.
CALL FOR FURTHER INFORMATION!

BUCKLE SHIELD™ Security Buckle available.
CALL FOR FURTHER INFORMATION!



IMPERVIOUS VINYL STRETCHER / COT STRAPS

METAL ROLLER BUCKLE

Item # 33071

Metal Roller Buckle, 7', 1 piece, each

Item # 33091

Metal Roller Buckle, 9', 1 piece, each

METAL DROP JAW BUCKLE

Item # 34091

Metal Drop Jaw Buckle, 9', 1 piece, each

Item # 34152

Metal Drop Jaw Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 34172

Metal Drop Jaw Buckle, w/ Loop Ends, 7', 2 piece, each

METAL DOUBLE D-RING

Item # 35091

Metal Double D-Ring Buckle, 9', 1 piece, each

Item # 350121

Metal Double D-Ring Buckle, 12', 1 piece, each

Item # 350151

Metal Double D-Ring Buckle, 15', 1 piece, each

PLASTIC LADDERLOC BUCKLE

Item # 36071

Plastic Ladderloc Buckle, 7', 1 piece, each

Item # 36091

Plastic Ladderloc Buckle, 9', 1 piece, each

PLASTIC SIDE RELEASE BUCKLE

Item # 37071

Plastic Side Release Buckle, 7', 1 piece, each

Item # 37091

Plastic Side Release Buckle, 9', 1 piece, each

Item # 37152

Plastic Side Release Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 37172

Plastic Side Release Buckle, w/ Loop Ends, 7', 2 piece, each

Item # 37252

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 37272

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 7', 2 piece, each

Item # 37352

Plastic Side Release Buckle, w/ Non-Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 37372

Plastic Side Release Buckle, w/ Non-Metal Swivel Speed Clip Ends, 7', 2 piece, each

Item # 37452

Plastic Side Release Buckle, w/ Big Brass Swivel Speed Clip Ends, 5', 2 piece, each

Item # 37472

Plastic Side Release Buckle, w/ Big Brass Swivel Speed Clip Ends, 7', 2 piece, each

PLASTIC DOUBLE ADJUST BUCKLE

Item # 38152

Plastic Double Adjust Buckle, w/ Loop Ends, 5', 2 piece, each

Item # 38172

Plastic Double Adjust Buckle, w/ Loop Ends, 7', 2 piece, each

Item # 38252

Plastic Double Adjust Buckle, w/ Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 38272

Plastic Double Adjust Buckle, w/ Metal Swivel Speed Clip Ends, 7', 2 piece, each

PLASTIC CAM BUCKLE

Item # 39071

Plastic Cam Buckle, 7', 1 piece, each

Item # 39091

Plastic Cam Buckle, 9', 1 piece, each

AVAILABLE IMPERVIOUS VINYL COLORS:

Black	Blue	Green
Orange	Yellow	



VINYL STICKY STRAPS

Item # 30051

Disposable Vinyl Sticky Straps, set of 3

Available Colors:

Black, Blue, Orange, Yellow and Green



Strap Carrying Case

The **Strap Carrying Case** is made with lined Heavy Duty Vinyl that is reinforced with 1" Nylon webbing. Velcro® closure. Belt loops on the back for easy carrying. Holds up to two straps.

Item #60004BK

Strap Carrying Case, each



SCREEN-PRINTING

Available in all webbing categories.

Call for further information.

614-444-2300 | www.dickmedical.com

CROSS HARNESS SYSTEM

NYLON

Item # 11170X

Nylon Cross Chest Harness, Shoulder & Chest

Item # 11170XSET

Nylon Cross Chest Harness Set, Shoulder, Chest, Leg & Ankle

***Set includes 1 each of item number 11152 & 1 each of item #11162**

Item # 11174XS

Nylon Cross Harness Shoulder Straps

Item # 11173XF

Nylon Cross Harness Female Strap, Chest X-Double Buckle

Item # 11173XM

Nylon Cross Harness Male Strap, Chest X-Buckle Tongue & Strap

BIOTHANE® G1

Item # 51170X

Biothane® G1 Cross Chest Harness, Shoulder & Chest

Item # 51170XSET

Biothane® G1 Cross Chest Harness Set, Shoulder, Chest, Leg & Ankle

***Set includes 1 each of item number 51152 & 1 each of item #51172**

Item # 51174XS

Biothane® G1 Cross Harness Shoulder Straps

Item # 51173XF

Biothane® G1 Cross Harness Female Strap, Chest X-Double Buckle

Item # 51173XM

Biothane® G1 Cross Harness Male Strap, Chest X-Buckle Tongue & Strap



CROSS HARNESS ONLY PICTURED



Hands Off™ Security Buckle available.
CALL FOR FURTHER INFORMATION!

BUCKLE SHIELD™ Security Buckle available.
CALL FOR FURTHER INFORMATION!



SCREEN-PRINTING

Available in all
webbing categories.

Call for further information.
614-444-2300
www.dickmedical.com

AVAILABLE BIOTHANE® COLORS:

Black	Blue	Golden Yellow
Neon Yellow	Orange	Red
Silver	Purple	

AVAILABLE NYLON COLORS:

Black	Blue	Green
Maroon	Olive Drab	Orange
Red	Yellow	

Meets SAE J3027 dynamic crash standard with Power-LOAD cot fastener system used with Power-PRO XT ambulance cot, Performance-LOAD cot fastener system used with Power-PRO XT ambulance cot, or Performance-LOAD cot fastener system used with Performance-PRO XT ambulance cot.

SHOULDER RESTRAINT HARNESS SYSTEMS

SHOULDER HARNESS RESTRAINT

The **Shoulder Harness Restraint System** is made up of two adjustable shoulder straps and one double adjusting torso strap. Metal Push Button buckle attaches into L-Links to secure on patient. Hands Off™ Security Buckle available. See various sizes below.

Item #11170C, 31170C, 51670C & 61170C fits Stryker Cot 6500 Series

- Available in **Nylon, Impervious, BioThane® G1**.

NYLON

Item # 11160

Nylon Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 5⅓' Torso Strap**

Item # 11153

Nylon Shoulder Harness Restraint System, Torso Strap **only**, each, **5⅓' Torso Strap**

Item # 11154

Nylon Shoulder Harness Restraint System, Shoulder Straps **only**, pair, **4½' Shoulder Straps**

Item #11170

Nylon Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 7⅓' Torso Strap**

Item #11170C

Nylon Shoulder Harness Restraint System, Complete Set, **6½' Shoulder Straps, 6' Torso Strap**

IMPERVIOUS

Item # 31160

Impervious Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 5⅓' Torso Strap**

Item # 31153

Impervious Shoulder Harness Restraint System, Torso Strap **only**, each, **5⅓' Torso Strap**

Item # 31154

Impervious Shoulder Harness Restraint System, Shoulder Straps **only**, pair, **4½' Shoulder Straps**

Item #31170

Impervious Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 7⅓' Torso Strap**

Item #31170C

Impervious Shoulder Harness Restraint System, Complete Set, **6½' Shoulder Straps, 6' Torso Strap**

BIO-THANE® G1

Item # 51660

BioThane® G1 Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 5⅓' Torso Strap**

Item # 51653

BioThane® G1 Shoulder Harness Restraint System, Torso Strap **only**, each, **5⅓' Torso Strap**

Item # 51654

BioThane® G1 Shoulder Harness Restraint System, Shoulder Straps **only**, pair, **4½' Shoulder Straps**

Item #51670

BioThane® G1 Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 7⅓' Torso Strap**

Item #51670C

BioThane® G1 Shoulder Harness Restraint System, Complete Set, **6½' Shoulder Straps, 6' Torso Strap**

SHOULDER HARNESS RESTRAINT WITH FIXED SHOULDER HOLDER

The Fixed Shoulder Holder is made with an 8" X 2" nylon strap sewn to both adjustable shoulder straps to prevent straps from falling off patients shoulders.

- Available in **Nylon, Impervious, BioThane® G1**.

Item # 11160F

Nylon Shoulder Harness System with Fixed Shoulder Holder, Complete Set

Item # 31160F

Impervious Shoulder Harness System with Fixed Shoulder Holder, Complete Set

Item # 51660F

BioThane® G1 Shoulder Harness System with Fixed Shoulder Holder, Complete Set

ADJUSTABLE SHOULDER HOLDER (ONLY)

Made with 1" black polypropylene webbing. Keeps shoulder straps in place.

Item # 24000BK

1" Black Polypropylene Shoulder Holder, each

NYLON



IMPERVIOUS



BIO-THANE®



BIOTHANE® G1

BioThane® G1 is Urethane coated polyester webbing that is dielectrically welded making it 100% impervious. BioThane® G1 straps are cleanable, waterproof and durable.

METAL PUSH BUTTON BUCKLE

Item # 51071

Metal Push Button Buckle, 7', 1 piece, each

Item # 51091

Metal Push Button Buckle, 9', 1 piece, each

Item # 51252

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 5', 2 piece, each

Item # 51272

Metal Push Button Buckle w/ Swivel Speed Clip Ends, 7', 2 piece, each

Item # 51652

Metal Push Button Buckle w/ Metal Loop Ends, 5', 2 piece, each

Item # 51672

Metal Push Button Buckle w/ Metal Loop Ends, 7', 2 piece, each

EXTENSION STRAPS

Item # 51021

Metal Push Button Buckle Extension Strap, 2', 1 piece, each

Item # 51031

Metal Push Button Buckle Extension Strap, 3', 1 piece, each

Item # 51041

Metal Push Button Buckle Extension Strap, 4', 1 piece, each

PLASTIC SIDE RELEASE BUCKLE

Item # 57071

Plastic Side Release Buckle, 7', 1 piece, each

Item # 57091

Plastic Side Release Buckle, 9', 1 piece, each

Item # 57252

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 5', 2 piece, each

Item # 57272

Plastic Side Release Buckle, w/ Metal Swivel Speed Clip Ends, 7', 2 piece, each

Item # 57652

Plastic Side Release Buckle, w/ Metal Loop Ends, 5', 2 piece, each

Item # 57672

Plastic Side Release Buckle, w/ Metal Loop Ends, 7', 2 piece, each

SHOULDER RESTRAINT HARNESS SYSTEMS

Item # 51660

BioThane® G1 Shoulder Harness Restraint System, Complete Set, 4½' Shoulder Straps, 5½' Torso Strap

Item # 51653

BioThane® G1 Shoulder Harness Restraint System, Torso Strap only, each, 5½' Torso Strap

Item # 51654

BioThane® G1 Shoulder Harness Restraint System, Shoulder Straps only, pair, 4½' Shoulder Straps

Item #51670

BioThane® G1 Shoulder Harness Restraint System, Complete Set, 4½' Shoulder Straps, 7½' Torso Strap

Item #51670C

BioThane® G1 Shoulder Harness Restraint System, Complete Set, 6½' Shoulder Straps, 6' Torso Strap



1B Metal Push Button Buckle	7B Plastic Side Release Buckle
2E Metal Swivel Speed Clip	6E Metal Loop End

AVAILABLE IN:

BLACK	BLUE
GOLDEN YELLOW	NEON YELLOW
ORANGE	RED
SILVER	



SCREEN-PRINTING
Available in all webbing categories.

Call for further information.
614-444-2300
www.dickmedical.com



Hands Off™ Security Buckle available.
CALL FOR FURTHER INFORMATION!

BUCKLE SHIELD™ Security Buckle available.
CALL FOR FURTHER INFORMATION!



G1 Webbing Break Strength: 1400 pounds

DISPOSABLE POLYPROPYLENE STRAPS

POLYPROPYLENE WEBBING BREAK STRENGTH: 1800 POUND BREAK STRENGTH

ECONOMY DISPOSABLE POLYPROPYLENE

Item # 47071

Plastic Side Release, 7', 1 piece, each

Item # 47091

Plastic Side Release, 9', 1 piece, each

Item # 47152

Plastic Side Release w/ Loop Ends, 5', 2 piece, each

Item # 47172

Plastic Side Release w/ Loop Ends, 7', 2 piece, each

Item # 47552

Plastic Side Release w/ Plastic Swivel Speed Clip Ends, 5', 2 piece, each

Item # 47572

Plastic Side Release w/ Plastic Swivel Speed Clip Ends, 7', 2 piece, each

Item # 48071

Plastic Double Adjust Buckle 7', 1 piece, each

Item # 48091

Plastic Double Adjust Buckle 9', 1 piece, each

Item # 48152

Plastic Double Adjust Buckle w/ Loop Ends, 5', 2 piece, each

Item # 48172

Plastic Double Adjust Buckle w/ Loop Ends, 7', 2 piece, each

Item #48552

Plastic Double Adjust Buckle w/ Plastic Swivel Speed Clip Ends, 5', 2 piece, each

Item # 48572

Plastic Double Adjust Buckle w/ Plastic Swivel Speed Clip Ends, 7', 2 piece, each

Item # 49071

Plastic Cam Buckle, 7', 1 piece, each

Item # 49091

Plastic Cam Buckle, 9', 1 piece, each

AVAILABLE POLYPROPYLENE COLORS:

BLACK	ORANGE
BLUE	RED
GREEN	YELLOW
NEON GREEN	WHITE



Plastic Side Release Buckle



Plastic Double Adjust Buckle



Loop End



Plastic Swivel Speed Clip

DISPOSABLE POLYPROPYLENE RED, WHITE & BLUE SETS

RED, WHITE & BLUE DISPOSABLE POLYPROPYLENE STRAP SETS

Item # 47071USA

Plastic Side Release, 7', 1 piece, set/3

Item # 47091USA

Plastic Side Release, 9', 1 piece, set/3

Item # 47152USA

Plastic Side Release w/ Loop Ends, 5', 2 piece, set/3

Item # 47172USA

Plastic Side Release w/ Loop Ends, 7', 2 piece, set/3

Item # 47552USA

Plastic Side Release w/ Plastic Swivel Speed Clip Ends, 5', 2 piece, set/3

Item # 47572USA

Plastic Side Release w/ Plastic Swivel Speed Clip Ends, 7', 2 piece, set/3

Item # 48071USA

Plastic Double Adjust Buckle 7', 1 piece, set/3

Item # 48091USA

Plastic Double Adjust Buckle 9', 1 piece, set/3

Item # 48152USA

Plastic Double Adjust Buckle w/ Loop Ends, 5', 2 piece, set/3

Item # 48172USA

Plastic Double Adjust Buckle w/ Loop Ends, 7', 2 piece, set/3

Item # 48552USA

Plastic Double Adjust Buckle w/ Plastic Swivel Speed Clip Ends, 5', 2 piece, set/3

Item # 48572USA

Plastic Double Adjust Buckle w/ Plastic Swivel Speed Clip Ends, 7', 2 piece, set/3

Item # 49071USA

Plastic Cam Buckle, 7', 1 piece, set/3

Item # 49091USA

Plastic Cam Buckle, 9', 1 piece, set/3

SETS OF

RED
WHITE
BLUE



Plastic Side Release Buckle



Plastic Cam Buckle



Plastic Double Adjust Buckle

SCREEN-PRINTING

Available in all webbing categories.

Call for further information.

614-444-2300 | www.dickmedical.com



BODY STRAP SYSTEMS

U BODY™ ECONOMY STRAP SYSTEM (BLACK)

The **Economy U Body™ Strap System** is made with 2" black polypropylene webbing. It has two (2) suspender style shoulder straps with four (4) body straps (chest, waist, leg and ankle). Hook and Loop closures allow for secure attachment.



Item # 20000BK

Economy Polypropylene U Body™ Strap, each

U BODY™ NYLON STRAP SYSTEM (MULTI-COLORED)

The **Multi-Colored U Body™ Strap System** is made with 2" nylon webbing. It has two (2) suspender style shoulder straps with five (5) color-coded body straps. Hook and Loop closures allow for secure attachment.




Item # 20100MC

Nylon Multi-Colored U Body™ Strap, each

Y BODY™ NYLON STRAP SYSTEM (MULTI-COLORED)

The **Multi-Colored Y Body™ Strap System** is made with 2" nylon webbing. It has two (2) Y style shoulder straps with four (4) body straps. Hook and Loop closures allow for secure attachment.



Item # 21000MC

Nylon Multi-Colored Y Body™ Strap, each

Item # 21000MC/CS

Nylon Multi-Colored Y Body™ Strap with Carrying Case, each

Y BODY™ NYLON STRAP SYSTEM (BLACK & RED)

The **Black & Red Y Body™ Strap System** is made with 2" nylon webbing. It has two (2) Y style shoulder straps with four (4) body straps. Hook and Loop closures allow for secure attachment.



Item # 21000BK/RD

Nylon Black/Red Y Body™ Strap, each

Item # 21000BK/RD/CS

Nylon Black/Red Y Body™ Strap with Carrying Case, each

Item # 21000BK

Nylon Black Y Body™ Strap, each (Not pictured)

Item # 21000BK/CS

Nylon Black Y Body™ Strap with Carrying Case, each (Not pictured)



DISPOSABLE BODY STRAP (DiBS™)

The **Disposable Body Strap System** is made with 2" vinyl webbing with industrial double sided sticky tape for a secure attachment.

Item # 33000

Disposable Body Strap (DiBS™), each

SCREEN-PRINTING

Available on **ALL** body straps.
Call for further information.

614-444-2300 | www.dickmedical.com



BODY STRAP SYSTEMS

Y BODY™ NYLON STRAP SYSTEM (BLUE & ORANGE)

The **Blue & Orange Y Body™ Strap System** is made with 2" nylon webbing. It has two (2) Y style shoulder straps with four (4) body straps. Hook and Loop closures allow for a secure attachment.

Item # 21000BL/OR

Nylon Blue/Orange Y Body™ Strap, each

Item # 21000BL/OR/CS

Nylon Blue/Orange Y Body™ Strap with Carrying Case, each

Item # 21000OR

Nylon Orange Y Body™ Strap, each (Not pictured)

Item # 21000OR/CS

Nylon Orange Y Body™ Strap with Carrying Case, each (Not pictured)



Y BODY™ POLYPROPYLENE STRAP SYSTEM (BLACK)


The **Black Y Body™ Strap System** is made with 2" polypropylene webbing. It has two (2) Y style shoulder straps with four (4) body straps. Hook and Loop closures allow for a secure attachment.

Item # 22000BK

Polypropylene Black Y Body™ Strap, each

Item # 22000BK/CS

Polypropylene Black Y Body™ Strap with Carrying Case, each



Y BODY™ POLYPROPYLENE STRAP SYSTEM (MULTI-COLORED)


The **Multi-Colored Y Body™ Strap System** is made with 2" polypropylene webbing. It has two (2) Y style shoulder straps with four (4) body straps. Hook and Loop closures allow for a secure attachment.

Item # 22000MC

Polypropylene Multi-Colored Y Body™ Strap, each

Item # 22000MC/CS

Polypropylene Multi-Colored Y Body™ Strap with Carrying Case, each



Y BODY™ POLYPROPYLENE STRAP SYSTEM (BLACK & RED)

The **Black & Red Y Body™ Strap System** is made with 2" polypropylene webbing. It has two (2) Y style shoulder straps with four (4) body straps. Hook and Loop closures allow for a secure attachment.

Item # 22000BK/RD

Polypropylene Black/Red Y Body™ Strap, each

Item # 22000BK/RD/CS

Polypropylene Black/Red Y Body™ Strap with Carrying Case, each



RED CORDURA BODY STRAP CARRYING CASE WITH ZIP CLOSURE

The Body Strap Carrying Case is an easy way to carry and store all Body Strap Systems. Made with heavy duty red Cordura with zip closure. *(fits all Body Strap styles)*

Screen Printing Available

Item #22000CS

Red Cordura Body Strap Carrying Case, each



DMS BRAND KITS

EMERGENCY DISPOSABLE OBSTETRICAL - STANDARD



Item #0112OB

Each kit contains:

- 1 Pair Non-Latex Exam Sterile Gloves
- 1 Each Disposable Scalpel, Sterile
- 1 Each OB Underpad, 17"x24"
- 2 Pack Gauze Sponge 4"x4", Sterile (2/pack)
- 2 Each Umbilical Clamp, Sterile
- 1 Each Placenta Bag 8"x15" with twist tie
- 1 Each OB Pad, Sterile
- 1 Each Receiving Blanket, 40"x48"

- 3 Each Disposable Non-Woven Towel 13"x18"
- 1 Each Bulb Syringe, Sterile
- 1 Each Plastic Apron
- 1 Each BZK Towelette
- 1 Each OB Towelette
- 1 Each Trash Bag
- 1 Each Apgar Scoring Chart

EMERGENCY DISPOSABLE OBSTETRICAL - WITH BABY BUNTING



Item #0113OB

Each kit contains:

- 1 Pair Non-Latex Exam Sterile Gloves
- 1 Each Disposable Scalpel, Sterile
- 1 Each OB Underpad, 17"x24"
- 2 Pack Gauze Sponge 4"x4", Sterile (2/pack)
- 2 Each Umbilical Clamp, Sterile
- 1 Each Placenta Bag 8"x15" with twist tie
- 1 Each OB Pad, Sterile
- 1 Each Receiving Blanket, 40"x48"
- 3 Each Disposable Non-Woven Towel 13"x18"

- 1 Each Bulb Syringe, Sterile
- 1 Each Plastic Apron
- 1 Each BZK Towelette
- 1 Each OB Towelette
- 1 Each Sterile Foil Baby Bunting
- 1 Each Trash Bag
- 1 Each Apgar Scoring Chart

EMERGENCY DISPOSABLE OBSTETRICAL - WITH BABY BLANKET & CAP



Item #0112OB-BC

Each kit contains:

- 1 Pair Sterile Non-Latex Exam Gloves
- 1 Each Disposable Scalpel
- 1 Each Receiving Blanket, 40" x 48"
- 2 Each Umbilical Clamps, Sterile
- 4 Each Disposable Non-Woven Towels, 13"x 18"
- 2 Each OB Towelettes
- 1 Each Bulb Syringe, Sterile
- 2 Each Alcohol Prep Pads

- 2 Each O.B. Pad, Sterile
- 1 Each Plastic bag for Placenta
- 2 Each Twist Ties
- 1 Each OB Underpad, 17" x 24"
- 3 Each Gauze Sponges, Sterile, Pack/2
- 1 Each Knit Baby Cap
- 1 Each Cloth Baby Blanket
- 1 Each Trash Bag
- 1 Each Apgar Scoring Chart

EMERGENCY DISPOSABLE OBSTETRICAL - WITH SAFETY GLASSES



Item #0113OB-NJ

Each kit contains:

- 4 Pair Non-Latex Exam Sterile Gloves
- 1 Each Scissors with Suture Kit, Sterile
- 6 Each Gauze Sponge 4"x4", Sterile (2/pack)
- 4 Each Umbilical Clamp, Sterile
- 2 Each Placenta Bag 8"x15" with twist ties
- 1 Each Receiving Blanket, 40"x48"
- 1 Each Safety Glasses
- 4 Each Disposable Non-Woven Towel 13"x18"
- 1 Each Bulb Syringe, Sterile

- 2 Each Alcohol Prep Pads
- 2 Each OB Towelette
- 1 Each Trash Bag
- 1 Each Apgar Scoring Chart

AMP Sterile Foil Baby Bunting



Item #9030SS

Each kit contains:

Reflective hooded foil blanket to conserve heat in newborn infants. (SINGLE USE)

Directions:

1. Open bunting and lay flat on work station
2. Place infant face up on bunting with head in hood
3. Remove adhesive protection strip
4. Wrap bunting around infant and secure with tape

Monitor infant temperature closely to prevent overheating.



DMS BRAND KITS AND ACCSEORIES

Biohazard Spill Kit

Item #6006

Each kit contains:

- | | |
|--|-------------------------------|
| 1 Each Disposable Gown with sleeves | 1 Each 4" Scoop / Scraper |
| 1 Each Disinfectant Cloth | 1 Each Pick Up Powder |
| 1 Each BZK Towelette | 3 Each Paper Towels |
| 1 Pair Latex Free, Powder Free Glove, Pair | 1 Each Antimicrobial Wipe |
| 1 Each Earloop face mask with eye shield | 1 Each E.Z. Zip Bag |
| 1 Each Biohazard Bag w/ Twist Tie | 1 Each Trash Bag w/ Twist Tie |



Biohazard Personal Protection Kit

Item #6005

Each kit contains:

- 1 Each Disposable Gown with Sleeves
- 1 Each Disinfectant Cloth
- 1 Each BZK Towelette
- 1 Pair Latex Free, Powder Free Gloves
- 1 Each Earloop face mask with eye shield
- 1 Each Biohazard Bag
- 1 Each 4" Scoop / Scraper
- 1 Each Pick Up Powder



Spill Clean Up Kit

Item #0119SK

Each kit contains:

- | | |
|---|------------------------------|
| 1 Pair Non-Latex Exam Gloves | 1 Each Biohazard Bag |
| 1 Each BZK Towelette | 1 Each Trash Bag |
| 1 Each Scoop and Scraper | 2 Each Twist Ties |
| 1 Each Red-Z Fluid Control Pick Up Powder | 1 Each 6"x12" Resealable Bag |
| 1 Each P.A.W.S. Antimicrobial Wipe | |
| 3 Each C-Fold Towels | |



Isolation Kit

Item #0117IK

Each Kit Contains:

- 1 Each Disposable Gown with Sleeves
- 1 Each Nonwoven Head Cover
- 1 Each Nonwoven Shoe Cover
- 1 Pair Latex Free, Powder Free Gloves
- 1 Each BZK Towelette
- 1 Each Biohazard Bag
- 1 Each Earloop Face Mask with Eye Shield



AMP Burn Sheet

Item #91039BS

DMS Burn Sheet

Sterile, tear-resistant, and individually sealed first aid covering that helps protect burn victims from contamination and additional trauma.

- Sterile
- Latex Free
- 60" X 90"



PPE

3-Ply FACE MASKS

Item #91019P3 - 3-Ply Face Mask, 50/box, Size: Universal, 50/BX

Item #91019P3L2 - Level 2 3-Ply Face Mask, Size: Universal, 50/BX

Item #91019P3CN - Children's 3-Ply Face Mask, Size: 5.7"X3.2", 50/BX

- Flat, easy on/easy off mask with standard pleats.
- Lightweight, comfortable and easy to breathe through.
- Bi-directional pleats for increased facial comfort and fit .

For use in environments with minimal fluid exposure. General patient care and other environments with minimum fluid exposure, isolation, IC practices, central processing, SPD, clinics, long term care facilities and home health care.

KN95 RESPIRATOR FACE MASK

Item #91019FM

50/box

Size: Universal (14.5 cm X 10.5 cm)

Standard: GB2626-2006

Type: ear-loop 4-ply

Size: 10.5 × 14.5 cm

Filter efficiency >95%

Shelf life: 3 years

Application: Anti-Virus Protection

FACE SHIELD - MEDICAL ISOLATION MASK

Regular - Item #91020FS

The face shield's wrap-around design provides the wearer with great peripheral vision while acting as a protective shield against non-hazardous liquid splash and light particles in a lab or controlled environment. The face shield is not ANSI Z87 rated for eye protection. You can wear the lightweight shield comfortably for long periods and it fits easily over prescription glasses or protective eyewear. For safe working conditions, the face shield is low distorting and 100% optically clear while the headband securely positions the device on the wearer. Clear, anti-fog protection with 180 degree coverage. Sold by EACH.

ISOLATION GOWNS

Item #91023CPE1M

CPE gown w/ thumb loops, AAMI Level 2 standard, Material: 100% polyethylene, One size fits most, Non-sterile, CASE/100

Item #91025PPPE1M

Isolation Gown Blue, LVL 1, PP+PE 40g, Over the head style, w/ thumb loops, Latex Free, Back Tie, CASE/100

Item #91025PPPE2M

AAMI Level 2 Gown, 30GSM, PP+PE Knitted Fabric Cuff, CASE/100

Item #91025PPPE3M

AAMI Level 3 Gown, 35GSM, PP+PE Knitted Fabric Cuff, CASE/100

Item #91025SMS4M

AAMI Level 4 Standard Isolation Gown, 45GSM, SMS Non-Woven with Lamination, Knitted Cuff, Sewn Seams, Size: XL, CASE/40

INFRARED THERMOMETERS

Item #91027IRT

Non-Contact Infrared Thermometer, Each

- Designed to test forehead temperature, which can adjust according to environmental temperature and forehead temperature.
- Adopts high-precise infrared sensor and has reliable performance.
- Saves 32 testing figures automatically.
- Large LCD display with backlit.
- Fahrenheit and centigrade available.
- Automatic switch-off.
- Batteries not included.



PPE

NITRILE EXAMINATION GLOVES

Item #91021NBLS

Nitrile Exam Gloves, Powder Free, Small, CASE/1000

Item #91021NBLM

Nitrile Exam Gloves, Powder Free, Medium, CASE/1000

Item #91021NBLL

Nitrile Exam Gloves, Powder Free, Large, CASE/1000

Item #91021NBLX

Nitrile Exam Gloves, Powder Free, X-Large, CASE/1000



VINYL EXAMINATION GLOVES

Item #91021VBLS

Vinyl Exam Gloves, Powder Free, Small, CASE/1000

Item #91021VBLM

Vinyl Exam Gloves, Powder Free, Medium, CASE/1000

Item #91021VBLL

Vinyl Exam Gloves, Powder Free, Large, CASE/1000

Item #91021VBLXL

Vinyl Exam Gloves, Powder Free, X-Large, CASE/1000



CAPS & COVERS

CAPS & COVERS

Item #91022BC4

Bouffant Caps, CASE/100

Item #91022SC - Not Pictured

Surgeon Cap w/Tie, BOX/500

Item #91024HN - Not Pictured

Hair Net, 10g PP Non-Woven, 21", CASE/2000



Item #91024BC

Beard Cover, CASE/2000

These disposable beard covers meet the FDA requirements for food handling and processing and are long enough to contain most beards in any sterile environment.



Item #91024SC

Disposable Shoe Covers, CASE/1800

These durable, non-slip covers are designed to keep dirt and contaminants out of your workplace, ensuring a clean and safe environment.




Item #91024FHC

Face and Head Cover, 10gsm, CASE/1000

Made from non-woven, breathable spun bound polypropylene fabric, these disposable, latex-free Balaclava hood covers are ideal for any food service or medical application.



BACKBOARDS



MASS CASUALTY SPINE BOARD

Item # 30003

The Mass Casualty spine board is a strong, yet economical choice, similar to our 10-993 spine board. This board is manufactured with mass casualty situations in mind, although it can be used in everyday emergency situations as well. This board comes with 10 pins for strapping and 21 large hand holes to make carrying easy. This board is not recommended for water usage. Load capacity is 600 lbs. Color: Yellow, Blue
Dimensions: 72.75" x 17.75" x 2.5" Weight: 18 lbs.
Available in Blue and Yellow

18" AB SPINE BOARDS

Item # 30002LG

The DMS AB spine board is an economical yet extremely strong board. The shell of the board is heavy duty and is x-ray translucent. The board comes with 21 handholds and 8 pins for your strapping needs. The spine board is waterproof and excellent for water rescue. The weight capacity of this board is 700lbs. It's made of durable polyethylene. In large orders, we can do custom colors. Colors: Lime Green, Orange, Red, Royal Blue, White, Yellow. QTY: 1 Product Dimensions: 72.75" x 17.75" x 2.5" Product Weight: 16.1 lbs.

Available in Lime Green, Orange, Red, Royal Blue, White and Yellow

PEDIATRIC SPINE BOARD (Not Pictured)

Item # 30004OR

The Pediatric Spine board is 100% X-Ray Translucent and is a must for camps and schools where there are children. We can't stress this enough, but a pediatric spine board is extremely important around pools and any area with many children because they normally get hurt while playing roughly. You can't put a small child on a 6-ft adult spine board properly.

Color: Orange

Product Dimensions: 53.5" x 14.75" x 2.5"

Product Weight: 10.45 lbs.

CHILD RESTRAINT SYSTEM

PEDIATRIC CHILD RESTRAINT SYSTEM

Item #57000BL

Pediatric Child Restraint System, each

- 3 point nylon strap system attaches to cot/stretcher
- 5 point nylon harness to secure patient
- Top cot strap measures 50"
- Plastic buckles
- Universal fitting - attaches to most cots and stretchers
- Fully adjustable straps
- Nylon carry case included
- Intended for patients 5-40 lbs.



CARRY POUCH™

CARRY POUCH™

Item #30090OR

Carry Pouch™, each

- Made with heavy orange vinyl.
- 1" black polypropylene webbing strap
- Rubberized and pre-molded handles to provide comfort
- Easy-to-clean
- Orange carrying case included



SOFT STRETCHERS

The **Scoop Up™ Stretcher** is designed to transport patients from confined areas. The patient is carried in an upright sitting position with 6 handholds. It measures 20" X 34" with 2" continuous and reinforced yellow or black nylon strapping. The Scoop Up™ can be easily stored. Made with Orange Vinyl makes it easy to clean .

Item #30070YL

Scoop Up™ Soft Stretcher, Yellow/Orange, each

Item #30070BK

Scoop Up™ Soft Stretcher, Black/Orange, each

***Comes with a handy carrying case.**



DMS **Soft Stretchers** are made with orange vinyl that is reinforced with black nylon strapping and vinyl handles.

- Rugged, yet comfortable vinyl handles added for easy gripping.
- Three sizes to choose from.

30" X 70" – **SMALL** Soft Stretcher - weight load - 1,000 pounds

40" X 71" – **MEDIUM** Soft Stretcher - weight load - 1,250 pounds

60" X 80" – **LARGE** Soft Stretcher - weight load - 1,500 pounds

Item # 30020OR

SMALL Soft Stretcher, each

Item # 30040OR

MEDIUM Soft Stretcher, each

Item # 30060OR

LARGE Soft Stretcher, each

Item #30001OR


Soft Stretcher Carrying Case, each



Carrying Case Included
and Sold Separately



LARGE



MEDIUM



SMALL

HIP & PELVIC STABILIZER SPLINTS



HIPPY HOLD™

The **HIPPY HOLD™** compression wrap helps stabilize the pelvic area. It will also help relieve pain in suspected fractures. Highly visible orange vinyl and color coded nylon straps with 2" adjustable plastic buckles provide ease of application.

Item #54000OR

Hippy Hold™, each

Item #54000ORM

Hippy Hold™, case/5



Carrying Case Included



HIP WRAP™

The **HIP WRAP™** is made with white corrugated plastic. This hip and pelvic stabilizer splint is a less painful way to stabilize and transport all hip and pelvic fractures. Adjustable straps are 2" polypropylene webbing with plastic side release buckles. Each Hip Wrap™ includes a 4" X 18" gray foam block. Comes with instruction sheet.

Item #53000

Hip Wrap™ each

Item #53000M

Hip Wrap™, case/12

POLYPROPYLENE, IMPERVIOUS VINYL & COTTON GAIT BELTS

COTTON GAIT BELTS

PLASTIC SIDE RELEASE BUCKLE

AVAILABLE COLORS:

Pastel Stripe (BPA), Natural (NAC), Light Pastel (LPA), Red & Blue Stripe (RBNA)

Item #67051

Cotton Gait Belt, Plastic Side Release Buckle, 60", each

Item #67061

Cotton Gait Belt, Plastic Side Release Buckle, 72", each

Item #67071

Cotton Gait Belt, Plastic Side Release Buckle, 84", each

METAL DROP JAW BUCKLE

AVAILABLE COLORS:

Pastel Stripe (BPA), Natural (NAC), Light Pastel (LPA), Red & Blue Stripe (RBNA)

Item #64051

Cotton Gait Belt, Metal Drop Jaw Buckle, 60", each

Item #64061

Cotton Gait Belt, Metal Drop Jaw Buckle, 72", each

Item #64071

Cotton Gait Belt, Metal Drop Jaw Buckle, 84", each



POLYPROPYLENE GAIT BELTS

PLASTIC SIDE RELEASE

AVAILABLE COLORS:

Black (BK), Blue (BL), Green (GR), Neon Green (NG)

Orange (OR), Red (RD), Yellow (YL), White (WH)

Item #27051

Polypropylene Gait Belt, Plastic Side Release Buckle, 60", each

Item #27061

Polypropylene Gait Belt, Plastic Side Release Buckle, 72", each

Item #27071

Polypropylene Gait Belt, Plastic Side Release Buckle, 84", each



METAL DROP JAW

AVAILABLE COLORS:

Black (BK), Blue (BL), Green (GR), Neon Green (NG)

Orange (OR), Red (RD), Yellow (YL), White (WH)

Item #24051

Polypropylene Gait Belt, Metal Drop Jaw Buckle, 60", each

Item #24061

Polypropylene Gait Belt, Metal Drop Jaw Buckle, 72", each

Item #24071

Polypropylene Gait Belt, Metal Drop Jaw Buckle, 84", each

IMPERVIOUS VINYL GAIT BELTS

PLASTIC SIDE RELEASE

AVAILABLE COLORS:

Black (BK), Blue (BL), Green (GR), Orange (OR), Yellow (YL)

Item #37051

Impervious Vinyl Gait Belt, Plastic Side Release Buckle, 60", each

Item #37061

Impervious Vinyl Gait Belt, Plastic Side Release Buckle, 72", each

Item #37071

Impervious Vinyl Gait Belt, Plastic Side Release Buckle, 84", each



METAL DROP JAW

AVAILABLE COLORS:

Black (BK), Blue (BL), Green (GR), Orange (OR), Yellow (YL)

Item #34051

Impervious Vinyl Gait Belt, Metal Drop Jaw Buckle, 60", each

Item #34061

Impervious Vinyl Gait Belt, Metal Drop Jaw Buckle, 72", each

Item #34071

Impervious Vinyl Gait Belt, Metal Drop Jaw Buckle, 84", each



BIOTHANE® G1 GAIT BELTS

BIOTHANE® GAIT BELTS

SPRING LOADED CAM

AVAILABLE COLORS:

**Black (BK), Blue (BL), Golden Yellow (GY), Neon Yellow (YL)
Orange (OR), Purple (PR), Red (RD), Silver (SL)**

Item #59051S

Biothane® G1 Gait Belt, Spring Loaded Cam Buckle, 60", each

Item #59061S

Biothane® G1 Gait Belt, Spring Loaded Cam Buckle, 72", each

Item #59071S

Biothane® G1 Gait Belt, Spring Loaded Cam Buckle, 84", each

PLASTIC SIDE RELEASE

AVAILABLE COLORS:

**Black (BK), Blue (BL), Golden Yellow (GY), Neon Yellow (YL)
Orange (OR), Purple (PR), Red (RD), Silver (SL)**

Item #57051

Biothane® G1 Gait Belt, Plastic Side Release Buckle, 60", each

Item #57061

Biothane® G1 Gait Belt, Plastic Side Release Buckle, 72", each

Item #57071

Biothane® G1 Gait Belt, Plastic Side Release Buckle, 84", each

METAL DROP JAW

AVAILABLE COLORS:

**Black (BK), Blue (BL), Golden Yellow (GY), Neon Yellow (YL)
Orange (OR), Purple (PR), Red (RD), Silver (SL)**

Item #54051

Biothane® G1 Gait Belt, Metal Drop Jaw Buckle, 60", each

Item #54061

Biothane® G1 Gait Belt, Metal Drop Jaw Buckle, 72", each

Item #54071

Biothane® G1 Gait Belt, Metal Drop Jaw Buckle, 84", each



BIOTHANE® G1

Made with **Biothane®** coated webbing, DMS gait belts are cleanable, waterproof and durable. They do not require machine washing and return to service 9X faster than traditional gait belts. With a antimicrobial coating, these gait belts resist bacteria, mold, and other microorganisms.

- Smooth strap easily wipes clean
- Available in plastic or metal buckle
- All edges are fully sealed and there are no seams
- Plastic buckle straps are MRI-safe



LIMB RESTRAINTS

NYLON WRIST and ANKLE RESTRAINTS - REUSABLE



Nylon Wrist Restraints are made of 2" nylon webbing with attached D-ring. Restraint wraps around patient's wrist and attaches to backboard or cot with Hook & Loop.

Item # 51000
Nylon Wrist Restraint, pair

Available Colors:

BLACK	BLUE	GREEN	MAROON
OLIVE DRAB	ORANGE	RED	YELLOW



Nylon Ankle Restraints are made of 2" nylon webbing. Restraint wraps around patient's ankles with Hook & Loop and attaches to backboard or cot with Hook & Loop.

Item # 52000
Nylon Ankle Restraint, pair

Available Colors:

BLACK	BLUE	GREEN	MAROON
OLIVE DRAB	ORANGE	RED	YELLOW

HEAVY DUTY NYLON LIMB RESTRAINTS - REUSABLE



Can be used for ankle and/or wrist application

HEAVY DUTY Nylon Limb Restraints are made of 2" nylon webbing. Restraint wraps around patient's ankles or wrists with Hook & Loop and attaches to backboard or cot with Hook & Loop.

Item # 52000
Nylon Limb Restraint, pair

Available Colors:

BLACK	BLUE	GREEN	MAROON
OLIVE DRAB	ORANGE	RED	YELLOW

WRIST & ANKLE - DISPOSABLE



NON-D-RING

Wrist & Ankle Restraints are made of soft 3/8" X 11" foam. 1" hook secures restraint on patient. Strap is 1" X 29" white polyester.

Item # 501100
Economy Limb Restraint, pair
Item # 501100M
Economy Limb Restraint, case/48 pair



D-RING

Wrist & Ankle Restraints with D-Rings are made of soft 3/8" X 11" foam. D-Rings secure restraint on patient. Strap is 1" X 51" white polyester.

Item # 501110
Economy Limb Restraint with D Rings, pair
Item # 501110M
Economy Limb Restraint with D Rings, case/ 48 pair



PLASTIC SIDE RELEASE

Wrist & Ankle Restraints with Plastic Side Release Buckles are made of soft 3/8" X 11" foam. Plastic Buckles secure restraint on patient. Strap is 1" X 51" white polyester.

Item # 501111
Economy Limb Restraint with Plastic Side Release Buckle, pair
Item # 501111M
Economy Limb Restraint with Plastic Side Release Buckle, case/ 48 pair

STAIR CHAIR REPLACEMENT STRAPS

Stair Chair replacement straps are 2" Polypropylene webbing with Plastic Side Release buckles. Set of 2 straps wraps around patient and chair to secure. Call for further information.

Item #47162BKGRSET

Economy Poly, Plastic SR Bkl, 6.5 ft., 2 pc., Black & Green, Set

Item #47162BKBKSET

Economy Poly, Plastic SR Bkl, 6.5 ft., 2 pc., Black X 2, Set

Item #47051

Economy Poly, Plastic SR Bkl, 5 ft., 1 pc., each



CHAIR NOT INCLUDED



Strap Carrying Case

Made with lined Heavy Duty Vinyl that is reinforced with 1" Nylon webbing. Velcro® closure. Belt loops on the back for easy carrying. Holds up to two straps.

Item #60004BK
Strap Carrying Case, each

Also available
in many different
colors.

BLACK	ORANGE
BLUE	RED
GREEN	YELLOW
MAROON	OLIVE DRAB
NEON GREEN	WHITE

EQUIPMENT HOLDING STRAP

Made of 1" black polypropylene webbing, 26" long, with 1" plastic side release buckle. Comes with two (2) 1" screw down brackets that secures to surfaces.

Equipment Holding Straps can be customized.

Item # 27132BK

Equipment Holding Strap, each



CINCH STRAP

Our Cinch Strap is made of 1½" Velcro hook & loop and a white loop buckle. Strap wraps around patients wrist and/or ankles and cinches with hook & loop Velcro®. Total strap measures 18" long.

Item # 49000BK

Black Cinch Strap, each

Item #49000OR

Orange Cinch Strap, each

AVAILABLE NYLON COLORS:

Black	Orange
-------	--------



LEND-A-HAND STRETCHER GRIPS

Ease transporting patients through confined spaces. Constructed of 2" Nylon webbing and a black, rubberized no-slip grip, pre-molded handle. Comes in sets of 4 each. Available in all standard Nylon colors.

Item # 17000

2" Nylon Webbing with Black PVC Plastic Handle, set/4

AVAILABLE COLORS:

BK	BL	GR	MA
OR	RD	YL	OD



HANDS OFF™ SECURITY BUCKLE

HANDS OFF™ SECURITY METAL PUSH BUTTON BUCKLE

Outfit your strap with the
NEW DMS Hands Off™ Security Buckle.
Fits all strap configurations,



HEAD IMMOBILIZERS

HEAD STAY™ HEAD IMMOBILIZER



The **Head Stay™** is an easy way to stabilize head or neck injuries. The Head Stay™ is made of water resistant cardboard. One blue adhesive head strap included. Each packaged flat for easy storage.

- 15" X 5" Base
- 15" head stay panels - Both sides of immobilizer are 6" when folded
- Can be easily adjusted to fit most patients
- Blue center foam pad for head comfort
- Highly visible red sticky pull tabs that secure firmly to board

Item # 50000
Head Stay™, each

ORANGE BLOCKS™



Using the **Orange Blocks™** is a fast and convenient way of immobilizing a patient's head. Head Pillows are sealed in orange vinyl bags making them water resistant. Each pillow is 4" X 5" X 10".

- They can be used with DMS blue adhesive head & chin straps and/or alignment pad for immobilization.

Head & Chin Straps



Alignment Pad



Item # 50200 Orange Blocks™, pair	Item #50200S Orange Blocks™ with Adhesive Head Straps, pair	Item #50200P Orange Blocks™ with Alignment Pad, pair	Item #50200SP Orange Block™ with Adhesive Head Straps & Alignment Pad, pair
Item #50209 Orange Blocks™, case/9	Item #50209S Orange Blocks™ with Adhesive Head Straps, case/9	Item #50209P Orange Blocks™ with Alignment Pad, case/9	Item #50209SP Orange Blocks™ with Adhesive Head Straps & Alignment Pad, case/9
Item #50230 Orange Blocks™, case/30	Item #50230S Orange Blocks™ with Adhesive Head Straps, case/24	Item #50230P Orange Blocks™ with Alignment Pad, case/30	Item #50230SP Orange Blocks™ with Adhesive Head Straps & Alignment Pad, case/24

STICKY WEDGE™ FOAM BLOCKS



The **Sticky Wedge™ Foam Blocks** are medium density blocks that stabilize patients with head and/or neck injuries. Each block is backed with adhesive sticky tape which can stick to any type of backboard.

- Squeeze & Peel
- No sticky residue when removed.

They can be used with DMS blue adhesive head & chin straps for immobilization.

Item # 50300
Sticky Wedge™ Foam Blocks, pair

Item # 50312
Sticky Wedge™ Foam Blocks, case/12 pair

Item # 50350
Sticky Wedge™ Foam Blocks, case/50 pair



HEAD IMMOBILIZERS

HEAD GUARD™ HEAD IMMOBILIZER

The **Head Guard™** is a reusable custom tailored head immobilizer that fits most spine boards. The side blocks (held in place with extra strong Hook & Loop that works even in water) are adjustable, shifting to accommodate almost any angle. Large ear holes make it easy for patient evaluation. The head and chin straps hold the head securely in place. Vinyl dipped foam blocks provide firm but resilient support and a clean hygienic surface. Easily attaches to spine boards with all Hook & Loop strapping.

DIMENSIONS:

- 15 ½" X 10" Base
- 6" H X 9½" W Tapered Pillows
- 39" Head & Chin Straps
- 25" attachment straps

Item # 50100

Head Guard™ Head Immobilizer, each

Item # 50106

Head Guard™ Head Immobilizer, case/6



Head Guard™ Pillows, Base and Head & Chin Straps ALSO SOLD SEPERATELY!

Item #50101

Head Guard™ Head Immobilizer Pillows, pair

Item #50102

Head Guard™ Head Immobilizer Base, each

Item #50103

Head Guard™ Head Immobilizer Head & Chin Straps, pair

AVAILABLE COLORS:

Black	Blue
Orange	Red



HARE TRACTION SPLINT REPLACEMENT STRAPS

ADULT & PEDIATRIC SIZES

Item #600230

(A) Ratchet Traction Device Strap W/ S Hook, Adult/Pediatric, each

Item #600221

(B) Traction Splint, Ischial 2 Piece Strap, Pediatric, each

Item #600220

(B) Traction Splint, Ischial 2 Piece Strap, Adult, each

Item #600211

(C) Traction Splint, Leg Straps, Pediatric, set

Item #600210

(C) Traction Splint, Leg Straps, Adult, set

Item #600201

(D) Traction Splint, Ankle Hitch, Pediatric, each

Item #600200

(D) Traction Splint, Ankle

Item #600001

Traction Splint Strap Replace Set, Pediatric

Set Includes:

- 1 ea. - Ratchet Traction Device Strap W/ S Hook
- 1 ea. - Traction Splint, Ischial 2 Piece Strap, Pediatric
- 4 ea. - Traction Splint, Leg Straps, Pediatric
- 1 ea. - Traction Splint, Ankle Hitch, Pediatric

Item #600002

Traction Splint Strap Replace Set, Adult

Set Includes:

- 1 ea. - Ratchet Traction Device Strap W/ S Hook
- 1 ea. - Traction Splint, Ischial 2 Piece Strap, Adult
- 4 ea. - Traction Splint, Leg Straps, Adult
- 1 ea. - Traction Splint, Ankle Hitch, Adult



PADDED BOARD SPLINTS

All splints are made with 1/2" X 3" wood, sealed in reusable orange vinyl bags. 1/2" foam on one side of splint adds comfort for patient.

Item # 60001

Splint Carrying Case, Velcro® Closure, each

Item #60002

Splint Carrying Case, Zipper Closure, each

Item #60100

Set of 2 each 15", 36", & 54" Padded Board Splints with Carrying Case

Item # 60000

Set of 2 each 15", 36" & 54", set

Item # 60015

3"X 15" Padded Board Splint

Item # 60015M

3"X 15" Padded Board Splint case/150

Item # 60036

3"X 36" Padded Board Splint

Item # 60036M

3"X 36" Padded Board Splint case/75

Item # 60054


3"X 54" Padded Board Splint

Item # 60054M

3"X 54" Padded Board Splint case/50



Velcro & Zippered
Carrying Cases



CARDBOARD SPLINTS

All of our **cardboard splints** are tapered for a more accurate fit. The 24" and 36" splints fold at the bottom for hand or foot placement. Use tape to secure.

PLAIN CARDBOARD SPLINTS

Item # 61012M

12" Economy Plain Cardboard Folding Splint, case/36

Item # 61018M

18" Economy Plain Cardboard Folding Splint, case/36

Item # 61024M

24" Economy Plain Cardboard Folding Splint, case/36

Item # 61036M

36" Economy Plain Cardboard Folding Splint, case/36

CENTER FOAM CARDBOARD SPLINTS

Item # 61012MCF

12" Economy Center Foam Cardboard Folding Splint, case/36

Item # 61018MCF

18" Economy Center Foam Cardboard Folding Splint, case/36

Item # 61024MCF

24" Economy Center Foam Cardboard Folding Splint, case/36

Item # 61036MCF

36" Economy Center Foam Cardboard Folding Splint, case/36

3 PIECE FOAM CARDBOARD SPLINTS

Item # 61012M3P

12" Economy 3 Piece Foam Cardboard Folding Splint, case/36

Item # 61018M3P

18" Economy 3 Piece Foam Cardboard Folding Splint, case/36

Item # 61024M3P

24" Economy 3 Piece Foam Cardboard Folding Splint, case/36

Item # 61036M3P

36" Economy 3 Piece Foam Cardboard Folding Splint, case/36

FULL FOAM CARDBOARD SPLINTS

Item # 61012MFF

12" Economy Full Foam Cardboard Folding Splint, case/36

Item # 61018MFF

18" Economy Full Foam Cardboard Folding Splint, case/36

Item # 61024MFF

24" Economy Full Foam Cardboard Folding Splint, case/36

Item # 61036MFF


36" Economy Full Foam Cardboard Folding Splint, case/36

24"

18"


PLAIN

12"




36"

CENTER FOAM



3 PIECE FOAM



FULL FOAM



PREMIUM & ECONOMY COTS

PREMIUM COT

Item #FC211

Military Style Folding Cot, w Bag, each

- 72"x25"
- Fabric: 600denier oxford,
- 330lb Capacity
- Comes with carry bag
- Standards: Latex Free
- Certification: N/A
- Reusable: Yes
- Application: Universal



ECONOMY COT

Item #EC121

Economy Cot with bag, each

- 71"x24"
- Fabric: 600denier oxford,
- 330lb Capacity
- Comes with carry bag
- Standards: Latex Free
- Certification: N/A
- Reusable: Yes
- Application: Universal



RESCUE CHAIR

RESCUE CHAIR

Item #RC100

DMS Rescue Chair, each

The DMS Rescue Chair is designed to serve as a wheeled chair, a stair chair, and a flat stretcher, and it allows you to quickly and easily switch to different positions depending for any situation. Features posts at the head-end of the stretcher as well as carrying handles and 2 stationary wheels at the foot-end. Measures 19.5" x 82" and can carry 350lbs

- Fixed head-end posts serve as support legs that are in position as soon as the chair is unfolded
- Fixed foot-end wheels allow easy maneuvering
- Folding foot-end handles provide excellent grip
- Only weighs 20lbs



CUSTOM SCREEN PRINT AVAILABLE.
CONTACT US FOR MORE INFORMATION!

614-444-2300

CARCINOGEN CONTAINMENT GEAR BAG

CARCINOGEN CONTAINMENT GEAR BAG




Item #GB-50015

Carcinogen Containment Gear Bag, each

- Reduces exposure to carcinogens and off-gassing from contaminated gear.
- Dielectrically welded seams provide a water tight seal.
- Easily decontaminated.
- Industrial-grade, highly-durable TPU laminated polyester.
- Heavy-gauge zippers that do not corrode.
- TPU-coated handles prevent fluid absorption and can be wiped clean.
- Corrosion-resistant high-grade steel hardware.
- Separate compartments to reduce cross-contamination risks.
- Red with Gold Maltese Cross.
- 28" x 15" x 14" (5,880 cu. in. of storage)
- Lightweight - only 3.3 lbs.

CARGO NET

CARGO NET



Custom made to fit your application.

DMS safety and cargo nets can be made to any size.

Available with metal seat belt buckles, heavy duty plastic buckles or track fittings. Please call for pricing.

Item #30150CN

Cargo Net, each

IV ARMBOARDS

DISPOSABLE IV ARMBOARDS

- Strong double wall corrugated core
- Comfortable foam padding
- Dielectrically sealed
- Washable satin finish
- Meets FDA requirements

Item # 62204

2"x 4" Disposable IV Arm board, each

Item # 62204M50

2"x 4" Disposable IV Arm board, case/ 50

Item # 62206

2"x 6" Disposable IV Arm board, each

Item # 62206M50

2"x 6" Disposable IV Arm board, case/50

Item # 62209

2"x 9" Disposable IV Arm board, each

Item # 62209M50

2"x 9" Disposable IV Arm board, case/50

Item # 62309

3"x 9" Disposable IV Arm board, each

Item # 62309M50

3"x 9" Disposable IV Arm board, case/50

Item # 62309M100

3"x 9" Disposable IV Arm board, case/100

Item # 62318

3"x 18" Disposable IV Arm board, each

Item # 62318M50

3"x 18" Disposable IV Arm board, case/50

REUABLE IV ARMBOARDS

- Strong plywood core
- Comfortable foam padding
- Dielectrically sealed
- Washable satin finish
- Meets FDA requirements

Item # 63204

2"x 4" Reusable IV Arm board, each

Item # 63204M50

2"x 4" Reusable IV Arm board, case/ 50

Item # 63206

2"x 6" Reusable IV Arm board, each

Item # 63206M50

2"x 6" Reusable IV Arm board, case/50

Item # 63209

2"x 9" Reusable IV Arm board, each

Item # 63209M50

2"x 9" Reusable IV Arm board, case/50

Item # 63309

3"x 9" Reusable IV Arm board, each

Item # 63309M50

3"x 9" Reusable IV Arm board, case/50

Item # 63309M100

3"x 9" Reusable IV Arm board, case/100

Item # 63318

3"x 18" Reusable IV Arm board, each

Item # 63318M50

3"x 18" Reusable IV Arm board, case/50





From stretcher/cot straps to head immobilizers, Dick Medical Supply has unparalleled quality as well as the fastest lead times, with most orders being shipped within 24-48 hours of receipt. We offer the lowest prices available with customer service that can't be beat!


DAVID DICK - ddick@dms4ems.com

GREG KNEELAND - gkneeland@dms4ems.com

KATRINA COPELAND - kcopeland@dms4ems.com



**DICK MEDICAL SUPPLY
1946 CAMARO AVENUE
COLUMBUS, OHIO 43207
614-444-2300 | WWW.DICKMEDICAL.COM**



DICK MEDICAL SUPPLY

1946 Camaro Avenue, Columbus, OH 43207

614-444-2300 ★ dickmedical.com